

Immunofluorescence of Semi-Thin Sections

- 1) Cut semi-thin sections and place in cluster in center of coated slides.
- 2) Get rid of Sucrose! --- very high affinity for fluorescence.
 - A. Outline Sucrose Droplet with Diamond Scraper.
 - B. Place Slides in Humid Chamber.
 - C. Cover Generously w/ plain PBS (over Sucrose droplet).
 - D. Wash for ~15 min.
- 3) Make Ammonium Chloride (weak base) 0.05M in PBS, for quenching (blocking) free aldehyde groups (26.5 mg/10 ml).
- 4) Shake PBS off of slide --- into garbage. Cover sections w/ Ammonium Chloride 0.05M. Quench for 15 min.
- 5) Weigh out 100 mg BSA in 10 ml PBS (sonicate briefly) = 1% BSA solution. Filter w/ Whatman #1 filter. Then take a 50 ml metric cylinder and fill it w/ 45 ml PBS and 5 ml 1% BSA = 0.1% BSA/PBS.
- 6) Shake Ammonium Chloride off of slides and cover w/ drop of Wash Buffer (0.1% BSA/PBS). Let remain for 10 min. Repeat 2 more times (total wash = 30 min).
- 7) Dilute Antibodies in 1% BSA. Need at least 50 ul antibody dilution per slide.
- 8) Shake off Washing Buffer. Wipe everything dry on slide except area where sections are. Add antibody. Two hour incubation.
- 9) Rinse off w/ 0.1% BSA/PBS w/ pipet (not directly on sections). Shake off. Cover w/ 0.1% BSA/PBS. Leave for ten minutes. Repeat two more times. (Total washing step = 30 min.)
- 10) If bridge is needed it is added here, diluted in 0.1M PBS/BSA. Incubate for one hour.
- 11) Rinse the same as in step 9.
- 12) Dilute Fluorescence Reagents in 1% BSA/PBS. Spin dilution for 5 min at speed 13 in Beckman Microfuge. Shake off buffer. Wipe slides clean except for area where sections lie. Add fluorescence. One hour incubation in light protected box.
- 13) Rinse same as in step #8, but last wash--- third time, just PBS.
- 14) Add coverslip:
 - A. Shake and Wipe off Ab (like in step 7).
 - B. Add drop of glycerol/PBS/PPD.
 - C. Add 22x40 Coverslips.
 - D. Seal with Nailpolish.
- 13) Observe.